

PRESTO A85t Process system


High cooling capacities enable extremely low temperatures down to -85 °C possible. The high heating capacity, particularly with the A85t and the W85t, provides even more flexibility in the application.

The highly dynamic temperature control systems PRESTO are designed for high-precision temperature control for a wide range of applications such as reactor vessels or material stress tests. Moreover, by using efficient components, the process systems can compensate exothermic and endothermic reactions exceptionally fast. Permanent internal monitoring and self-lubricating pumps ensure a long life-time. In addition, numerous interfaces offer many remote control possibilities across networks or for integration into higher-level control systems.



Product features

- Cooling capacity up to 2.8 kW
- External Pt100 sensor connection
- · Connections for alarm-output, Pt100 external sensor and standby-input
- Analog connections, RS485, Profibus DP (accessory)
- Built-in 5.7" industrial color touchscreen
- Second external Pt100 sensor connection (accessory)
- Alarm output
- Heating capacity up to 15 kW
- Temperature stability ±0.05 °C ... ±0.1 °C
- Pump pressure up to 3 bar, max. flow rate 80 l/min
- removable ventilation grid







Performance values

400V/3PNPE/50Hz (Plug 32A CEE)	
Heating capacity kW	15
Viscosity max. cSt	50
Pump capacity flow rate I/min	0 80
Pump capacity flow pressure bar	0.1 3
Power consumption A	31

Order No.				9420852.07.T			
Cooling capacity (Ethanol)							
°C	20	0	-20	-30	-40	-60	-80
kW ¹	2.5	2.4	2.4	2.4	2.4	2.2	0.4
Refrigerant stage 1			Refrigerant stage 2				
Refrigerant		R507		Refrigerant		R23	
Filling weight g		1600		Filling weight g	3	680	
Global Warming Potential for 39 R507		3985		Global Warming Potential for R23		14800	
Carbon dioxide	equivalent t	6.376		Carbon dioxide	e equivalent t	10.064	

¹ Performance specifications measured in accordance with DIN 12876. Cooling capacities up to 20 °C measured with ethanol; over 20 °C with thermal oil unless otherwise specified. Performance specifications apply at an ambient temperature of 20 °C. Performance values may differ with other bath fluids.

Technical data

Available voltage versions		Cooling		
Order No.	Order No. 9 420 852.T		Cooling of compressor	2-stage Air
Available voltage versions:				
9420852.07.T 400V/3PNPE/50Hz (Plug 32A CEE) (R507)				
9420852.16.T 208-230V/3PPE/60Hz (Without Plug) (R507)				
9420852.06.T	230V/3PPE/50Hz	: (Plug 63A CEE) (R507)		
Other			Electronics	
Sound pressure level db	A	69	Interfaces	Alarm output, Ethernet,
Classification IP Code		Classification III (FL)		Modbus, Profibus optional, REG/EPROG
		IP 20	optional, RS232, RS485 optional, SD	
Pump typeCentrifugal PumpPump type Magnetically coupled1				memory card, Standby- Input optional, USB
			External pt100 sensor connection	integrated
			2nd external Pt100 sensor connection	accessory
			Integrated programmer	8x60 steps
			Temperature control	ICC
		Absolute temperature calibration	3 Point Calibration	
		Temperature display	5.7" TFT Touchscreen	
			Temperature setting	Touchscreen
Dimensions and vo	lumes		Temperature values	

Julaba

Internal usable expansion volume I	7	Setting the resolution of the temperature	0.01	
Minimal process volume l	9.5	display °C		
Active heat exchanger volume I	5	Working temperature range °C	-85 +250	
Weight kg	365	Temperature stability °C	±0.05 ±0.1	
5 5		Ambient temperature °C	+5 +40	
Dimensions cm (W × L × H)	61 x 108 x 125	Setting the resolution of the temperature	0.01	
Pump connections	M30x1.5 male	display °C	0.01	

All Benefits



Touch display. Perfect operation.

With the touch display, the user always has an overview of all values and functions. The intuitive and multilingual menu structure enables perfect control.



100 % Cooling capacity

'Active Cooling Control' for cooling available throughout the entire working temperature range, fast cool-down even at higher temperatures



Full control

'Temperature Control Features', for individual optimization, access to all important control parameters, additional settings for band limit, limits, co-speedfactor etc.



Highest measuring accuracy

'Absolute Temperature Calibration' for manual compensation of a temperature difference, 3-point calibration



Many interfaces.

Straight-forward remote control, data management, and integration into process structures. USB, Ethernet, RS232, SD card, and alarm off are permanently integrated. Further interfaces available as accessories.



Continuous operation up to +40 °C Robust temperature control instrument, continuous operation even at ambient temperatures of up to +40 °C



Duplicate safety

Adjustable high temperature cut-off for internal tank and for integrated expansion vessel



Convenience for several users

Administrator level for customizing instrument settings, user levels with limited permissions for fast and safe defined access, password protection, all levels adjustable



Intelligent temperature control. Intelligent cascade control - automatic and selfoptimizing adaptation of the PID control parameters with external stability of +/- 0.05 °C.



Control of the external application

External Pt100 sensor connection for precise measurement and control directly in the external application



Intelligent pump system

Reliable and consistent pump capacity, electronically adjustable pump stages or pressure value, automatic adjustment of pump capacity to viscosity



ax

m

Space-saving footprint

All connections as well supply and exhaust air are located at the front or rear, no venting grids on the sides, units can be placed close to each other or the application

Maximum safety.

Classification III according to DIN12876-1 enables safe operation, even with flammable fluids. Automatic switch-off in the event of high temperature or low liquid level.



For flammable bath fluid Classification III (FL) according to DIN 12876-1

3





Quick support

If an error occurs, the integrated Black-Box function permits fast diagnosis by the JULABO service team



Green technology. Development consistently applied environmentally

friendly materials and technologies.



Quick start.

Individual JULABO consultation and comprehensive manuals at your disposal.



100% Checked.

100% testing. 100% quality. Each JULABO Circulator undergoes thorough quality testing before leaving the factory.



8

JULABO. Quality.

Highest standards of quality for a long product life.

Satisfied customers.

11 subsidiaries and more than 100 partners worldwide guarantee fast and qualified JULABO support.



Services 24/7.

Around the clock availability. You can find suitable accessories, data sheets, manuals, case studies, and more at www.julabo.com.