

MAGIO MS-310F Refrigerated / heating circulator

As with all circulators from the MAGIO range, the refrigerated circulators stand out thanks to their premium quality, high performance and intuitive operation. The devices offer extra strong pressure and suction pumps, thus fulfilling the highest demands for temperature control of external applications. Whether in basic research, material testing or technical systems – the MAGIO refrigerated circulators offer high-tech solutions for high customer requirements.

With natural refrigerant

The MAGIO MS-310F with natural refrigerant. Order No.. 9 032 713.N1



Product features

- Ideal for demanding external applications
- Simple control of complex applications
- Continuously adjustable, extremely powerful pressure / suction pump
- Flow rate 16 ... 31 I / min, pressure 0.24 ... 0.92 bar, suction 0.03 ... 0.4 bar
- · Large, high-resolution TFT touch display with multilingual user interface
- · Parts being in contact with the medium made of stainless steel
- Integrated programmer
- Integrated external Pt100 connection
- USB connection
- RS232 interface for online communication
- Ethernet
- analog interfaces (accessory)
- Class III (FL) according to DIN 12876-1
- Connections for solenoid valve
- Integrated pump connection M16×1

Cool-down time Heat-up time °C 200 0 150 -20 100 310F -40 50 20 60 min Medium: Thermal Medium: Ethanol





Performance values

230V/50Hz (Schuko Plug - CEE 7/4 Plug Type F)		
Heating capacity kW	2	
Viscosity max. cSt	70	
Pump capacity flow rate I/min	16 31	
Pump capacity flow pressure bar	0.24 0.92	
Maximum suction bar	-0.030.4	
Power consumption A	14	

Refrigerant variants

Order No.	9032713.N1.33				
Cooling capacity 1 (Ethanol)				
°C	20	0	-10	-20	-30
kW ¹	0.33	0.28	0.23	0.13	0.05
Cooling capacity 2 (Ethanol)				
°C	20	0	-10	-20	-30
kW ¹	0.26	0.21	0.17	0.1	0.01

Note about natural refrigerants:

Temperature control units using natural refrigerants are often subject to regulatory requirements regarding the installation site, operation, transport or disposal of the units. If you have any questions, we will be happy to advise you.

Refrigerant stage 1	
Refrigerant	R290
Filling weight g	30
Global Warming Potential for R290	3
Carbon dioxide equivalent t	9.0E-5

Order No.	9032713.S1.33				
Cooling capacity 1 (I	Ethanol)				
°C	20	0	-10	-20	-30
kW ¹	0.33	0.28	0.23	0.13	0.05
Cooling capacity 2 (I	Ethanol)				
°C	20	0	-10	-20	-30
kW ¹	0.26	0.21	0.17	0.1	0.01
Refrigerant stage 1					

Refrigerant stage 1	
Refrigerant	R449A
Filling weight g	63

Global Warming Potential for R449A

1397

Carbon dioxide equivalent t 0.088011

Technical data

¹ Performance specifications measured in accordance with DIN 12876. Cooling capacities up to 20 °C measured with ethanol; over 20 °C with thermal oil unless otherwise specified. Performance specifications apply at an ambient temperature of 20 °C. Performance values may differ with other bath fluids. Cooling capacity 1 = capacity at minimum pump level, cooling capacity 2 = capacity at maximum pump level



Available voltage ver	SIONS
Order No.	9 032 713
Available voltage versions	:
9032713.N1.01	100V/50-60Hz (Nema N5-15 Plug) (R290)
9032713.S1.01	100V/50-60Hz (Nema N5-15 Plug) (R449A)
9032713.N1.02	115V/60Hz (Nema N5-15 Plug) (R290)
9032713.S1.02	115V/60Hz (Nema N5-15 Plug) (R449A)
9032713.N1.33	200-230V/50-60Hz (Schuko Plug - CEE 7/4 Plug Type F) (R290)
9032713.S1.33	200-230V/50-60Hz (Schuko Plug - CEE 7/4 Plug Type F) (R449A)
9032713.N1.33.chn	200-230V/50-60Hz (CN Plug) (R290)
9032713.S1.33.chn	200-230V/50-60Hz (CN Plug) (R449A)
9032713.N1.04	200-230V/50-60Hz (UK Plug Type BS1363A) (R290)
9032713.S1.04	200-230V/50-60Hz (UK Plug Type BS1363A) (R449A)
9032713.N1.05	200-230V/50-60Hz (CH Plug Type SEV 1011) (R290)
9032713.S1.05	200-230V/50-60Hz (CH Plug Type SEV 1011) (R449A)

Cooling	
Cooling of compressor	1-stage Air

Bath		Other	
Bath tank	Stainless steel	Classification	Classification III (FL)
Bath cover	integrated	IP Code	IP 20
Usable bath opening cm (W x L / D) 13 x 15 / 15	13 x 15 / 15	Pump function	Pressure Suction Pump
		Pump type	Immersion Pump
		User Interface Language	Chinese, Czech, Dutch, English, French, German, Italian, Japanese, Korean, Polish, Portuguese, Russian, Spanish, Turkish

	Dimensions and volumes	
Ethernet, Modbus,	Weight kg	29
RS232, RS485, Stakei, USB	Dimensions cm (W \times L \times H)	23 x 40 x 65
integrated	Filling volume I	3 4
8x60 steps	Pump connections	M16x1 male
ICC		
10 Point Calibration		
7" TFT Touchscreen		
Touchscreen		
00:00 99:59		
i 8	RS232, RS485, Stakei, JSB Integrated BX60 steps CC IO Point Calibration TTTT Touchscreen	Ethernet, Modbus, RS232, RS485, Stakei, JSB Dimensions cm (W × L × H) Filling volume I Pump connections CC TO Point Calibration T'' TFT Touchscreen Fouchscreen

Temperature values		Incl
Setting the resolution of the temperature display °C	0.01	2 Bar M16
Working temperature range °C	-30 +200	

Included in delivery

2 Barbed fittings for tubing 8 and 12 mm ID. (Pump connections M16x1 male)



Temperature stability °C	±0.01
Ambient temperature °C	+5 +40
Setting the resolution of the temperature display °C	0.01

All Benefits



Intelligent temperature control.

Intelligent cascade control - automatic and selfoptimizing adaptation of the PID control parameters with external stability of +/- 0.05 °C.



Many interfaces.

Straight-forward remote control, data management, and integration into process structures. USB, Ethernet, RS232, SD card, and alarm off are permanently integrated. Further interfaces available as accessories.



Touch display. Perfect operation.

With the touch display, the user always has an overview of all values and functions. The intuitive and multilingual menu structure enables perfect



Maximum safety.

Classification III according to DIN12876-1 enables safe operation, even with flammable fluids. Automatic switch-off in the event of high temperature or low liquid level.



Space saving. Free up space.

Place your JULABO Circulator right next to an application, another unit, or wall. That saves space. This is made possible by eliminating vents and connections on the sides.



Multi-lingual.

Operation in multiple languages.



Programmer. Integrated.

The integrated internal programmer makes it possible to automatically run temperature time profiles.



Analog I/O.

Analog interfaces for integration into process control systems (optional).



Temperature. Under control.

External Pt100 sensor connection for precise measurement and control directly in the external application.



Fill level. Monitored.

Fill level indicator on the display for heat-transfer liquid.



Process stability.

Early warning - visual and acoustic - of critical states increases process stability.



Process. Under control.

Full regulation of the dynamics control, access to all important control parameters for individual process optimization.



Stable. Mobile.

Rubber feet keep JULABO Circulators standing firm. Larger and more powerful units also have integrated rollers for easy handling.



Energy-saving.

The high-quality insulation of all relevant components saves energy.



Everything made of stainless steel.

Quality and material compatibility at the highest level. All parts in contact with the medium are entirely made of stainless steel.



Most powerful pump.

The integrated pressure/suction pump with performance values of 0.9 bar and -0.4 bar is the most powerful in its class and continuously adjustable.





Wide range.

Refrigerated and heating circulator in various combinations, circulator in various sizes.

Maximum flexibility through a large selection of accessories.



Condensation protection.

Superb design solution. Integrated ventilation directs air over the bath lid and minimizes condensation.



Connection. Easy.

Inclined pump connections (M16×1) facilitate the connection of applications. Each unit includes 2 barbed fittings of 8/12 mm diameter each.



Highest measuring accuracy

'Absolute Temperature Calibration' for manual compensation of a temperature difference, 10point calibration



100% Checked.

100% testing. 100% quality. Each JULABO Circulator undergoes thorough quality testing before leaving the factory.



Green technology.

Development consistently applied environmentally friendly materials and technologies.



Condensation protection.

Superb design solution. Integrated ventilation directs air over the bath lid and minimizes condensation.



Handle with ease.

Makes day-to-day work easy. Comfortably move your JULABO Circulator around by using the ergonomic handles (front and rear).



JULABO. Quality.

Highest standards of quality for a long product life.



Internal. External.

The pump is controlled via a lever located directly below the display. Easily change between internal and external circulation.



More bath.

Designed for more comfort. Thanks to the recessed cooling coil, the internal bath provides more space.



Quick start.

Individual JULABO consultation and comprehensive manuals at your disposal.



Satisfied customers.

11 subsidiaries and more than 100 partners worldwide guarantee fast and qualified JULABO support.



Solid

Minimized energy loss through high-quality insulation.



Services 24/7.

Around the clock availability. You can find suitable accessories, data sheets, manuals, case studies, and more at www.julabo.com.



Stable.

Rubber feet allow for a secured footing of your CORIO to prevent damage to your laboratory equipment.



Quick support

If an error occurs, the integrated Black-Box function permits fast diagnosis by the JULABO service team



Tidy.

The special drain tap for easy draining of bath fluids without tools.



Touching permitted.

Optimum safety. The ergonomic plastic handle protects your fingers from hot surfaces.