

# PRESTO A30 Process system

The PRESTO A30 is JULABO's smallest process system and combines high performance and a very compact design with all the advantages of the PRESTO series for a working temperature range from -30 °C to +250 °C.

The highly dynamic temperature control systems PRESTO are designed for high-precision temperature control for a wide range of applications such as reactor vessels or material stress tests. Moreover, by using efficient components, the process systems can compensate exothermic and endothermic reactions exceptionally fast. Permanent internal monitoring and self-lubricating pumps ensure a long life-time. In addition, numerous interfaces offer many remote control possibilities across networks or for integration into higher-level control systems.



## **Product features**

- Heating capacity up to 2.7 kW
- Temperature stability ±0.01 °C ... ±0.05 °C
- Alarm output
- Built-in 5.7" industrial color touchscreen
- External Pt100 sensor connection
- Second external Pt100 sensor connection (accessory)
- Connections for USB, Ethernet, RS232, and Alarm output
- Analog connections, RS485, Profibus DP (accessory)
- removable ventilation grid
- RS232
- RS485 (Accessory)
- Profibus DP (Accessory)

## Pump capacity



## Heat-up time



Medium: Thermal HL

## Cool-down time



## Medium: Thermal HL



## Performance values

| 230V/50Hz (Schuko Plug - CEE 7/4 Plug Type F) |     |  |
|---|-----|--|
| Heating capacity kW                           | 2.7 |  |
| Viscosity max. cSt                            | 50  |  |
| Pump capacity flow rate I/min                 | 25  |  |
| Pump capacity flow pressure bar               | 0.5 |  |
| Power consumption A                           | 15  |  |

| Order No.                          | 9420300.03 |     |     |     |      |
|------------------------------------|------------|-----|-----|-----|------|
| Cooling capacity 1 (Ethanol)       |            |     |     |     |      |
| °C                                 | 20         | 0   | -10 | -20 | -30  |
| kW <sup>1</sup>                    | 0.5        | 0.4 | 0.3 | 0.2 | 0.02 |
| Refrigerant stage 1                |            |     |     |     |      |
| Refrigerant                        | R452A      |     |     |     |      |
| Filling weight g                   | 340        |     |     |     |      |
| Global Warming Potential for R452A | 2140       |     |     |     |      |
| Carbon dioxide equivalent t        | 0.7276     |     |     |     |      |

<sup>1</sup> Performance specifications measured in accordance with DIN 12876. Cooling capacities up to 20 °C measured with ethanol; over 20 °C with thermal oil unless otherwise specified. Performance specifications apply at an ambient temperature of 20 °C. Performance values may differ with other bath fluids.

## **Technical data**

| Available voltage versions       |  | Cooling                              |  |                      |
|----------------------------------|--|--------------------------------------|--|----------------------|
| Order No.                        | 9 420 300                                |                                      | Cooling of compressor  | Air                  |
| Available voltage versions:      |  |                                      |  |                      |
| 9420300.03                       | 200-230V/50Hz (Sc<br>Plug Type F) (R452/ | huko Plug - CEE 7/4<br>A)            |  |                      |
| 9420300.04                       | 230V/50Hz (UK Plu<br>(R452A)             | g Type BS1363A)                      |  |                      |
| 9420300.05                       | 230V/50Hz (CH Plu<br>(R452A)             | g Type SEV 1011)                     |  |                      |
| 9420300.14                       | 208V/60Hz (Nema                          | N6-20 Plug) (R449A)                  |  |                      |
| Other                            |  |                                      | Electronics  |                      |
| Sound pressure level d           | lbA                                      | 54                                   | Interfaces Alarm output, E<br>Modbus, Profit<br>optional, REG/<br>optional, RS23<br>RS485 optional<br>memory card, S |                      |
| Classification                   |  | Classification III (FL)              |  |                      |
| IP Code                          |  | IP 21                                |  |                      |
| Pump type                        |  | Centrifugal Pump                     |  |                      |
| Pump type Magnetically coupled 1 |  | 1                                    |  | Input optional, USB  |
|                                  |  |                                      | External pt100 sensor connection   | integrated           |
|                                  |  | 2nd external Pt100 sensor connection | not available  |                      |
|                                  |  | Integrated programmer                | 8x60 steps   |                      |
|                                  |  |                                      | Temperature control  | ICC                  |
|                                  |  |                                      | Absolute temperature calibration   | 3 Point Calibration  |
|                                  |  |                                      | Temperature display  | 5.7" TFT Touchscreen |

Temperature setting

Touchscreen

# lulabo

| Dimensions and volumes                  |              |
|---|--------------|
| Internal usable expansion volume I      | 1.5          |
| Minimal process volume l                | 2.4          |
| Active heat exchanger volume I          | 1.4          |
| Weight kg                               | 62           |
| Dimensions cm ( $W \times L \times H$ ) | 25 x 59 x 62 |
| Pump connections                        | M24x1.5 male |

| Temperature values                                   |             |
|--|-------------|
| Setting the resolution of the temperature display °C | 0.01        |
| Working temperature range °C                         | -30 +250    |
| Temperature stability °C                             | ±0.01 ±0.05 |
| Ambient temperature °C                               | +5 +40      |
| Setting the resolution of the temperature display °C | 0.01        |

# All Benefits



## Convenience for several users

Administrator level for customizing instrument settings, user levels with limited permissions for fast and safe defined access, password protection, all levels adjustable



## Intelligent temperature control.

Intelligent cascade control - automatic and selfoptimizing adaptation of the PID control parameters with external stability of +/- 0.05 °C.



## Control of the external application

External Pt100 sensor connection for precise measurement and control directly in the external application



## Intelligent pump system

Reliable and consistent pump capacity, electronically adjustable pump stages or pressure value, automatic adjustment of pump capacity to viscosity



## Space-saving footprint

All connections as well supply and exhaust air are located at the front or rear, no venting grids on the sides, units can be placed close to each other or the application



## Maximum safety.

Classification III according to DIN12876-1 enables safe operation, even with flammable fluids. Automatic switch-off in the event of high temperature or low liquid level.



For flammable bath fluid Classification III (FL) according to DIN 12876-1



## 100 % Cooling capacity

'Active Cooling Control' for cooling available throughout the entire working temperature range, fast cool-down even at higher temperatures



## Full control

'Temperature Control Features', for individual optimization, access to all important control parameters, additional settings for band limit, limits, co-speedfactor etc.



## Highest measuring accuracy

'Absolute Temperature Calibration' for manual compensation of a temperature difference, 3-point calibration



#### Straight-forward remote control, data management, and integration into process structures. USB, Ethernet, RS232, SD card, and alarm off are permanently integrated. Further interfaces available as accessories.



Continuous operation up to +40 °C Robust temperature control instrument, continuous operation even at ambient temperatures of up to +40 °C



**Duplicate safety** Adjustable high temperature cut-off for internal tank and for integrated expansion vessel

# Quick support

If an error occurs, the integrated Black-Box function permits fast diagnosis by the JULABO service team



3





100% Checked.

100% testing. 100% quality. Each JULABO Circulator undergoes thorough quality testing before leaving the factory.



JULABO. Quality. Highest standards of quality for a long product life.



# Quick start.

Individual JULABO consultation and comprehensive manuals at your disposal.



## Services 24/7.

Around the clock availability. You can find suitable accessories, data sheets, manuals, case studies, and more at www.julabo.com.



# Green technology.

Development consistently applied environmentally friendly materials and technologies.



## Touch display. Perfect operation.

With the touch display, the user always has an overview of all values and functions. The intuitive and multilingual menu structure enables perfect control.

#### Satisfied customers.

11 subsidiaries and more than 100 partners worldwide guarantee fast and qualified JULABO support.